

ABUTMENT SCREW

TORQUE AND DRIVER GUIDE

2018

MANUFACTURER

TORQUE N-CM

DRIVER

| MANUFACTURER | TORQUE N-CM | | | | | | | | DRIVER |
|---|-------------|----|----|----|----|----|----|----|-------------------------------------|
| | 15 | 20 | 24 | 25 | 30 | 32 | 35 | 45 | |
| Astra 3.0 | X | | | | | | | | 0.050"-1.25mm hex |
| Astra 3.5-4.0 | | X | | | | | | | 0.050"-1.25mm hex |
| Astra 4.5-5.0 | | | | X | | | | | 0.050"-1.25mm hex |
| Astra EV 3.0, 3.6, 4.2, 4.8, 5.4 | | | | X | | | | | 0.050"-1.25mm hex |
| BioHorizons® External, Internal, Internal Tapered | | | | | X | | | | 0.050"-1.25mm hex |
| Biomet 3i™ External Hex | | | | | | | X | | 0.048" hex |
| Biomet 3i™ Certain® | | X | | | | | | | 0.048" hex |
| BlueSkyBio One Stage | | | | | X | | | | Star & 0.048" hex |
| BlueSkyBio Trilobe | | | | | X | | | | Unigrip & 0.048" hex |
| BlueSkyBio Internal Hex & Molar | | | | | X | | | | 0.050" hex |
| BlueSkyBio Conus 12 & Three | | | | | X | | | | 0.050" hex |
| BlueSkyBio Max | | | | | X | | | | 0.048" hex |
| BlueSkyBio Quatro | | | | | X | | | | 0.048" hex |
| Brånemark (Nobel Biocare) | | | | | | | X | | Unigrip |
| Camlog™ & Conelog® | | X | | | | | | | 0.050"-1.25mm hex |
| Dentium Super Line & Implantium | | | | | X | | | | 1.28mm hex |
| Dentsply Friadent/Frialit & Xive | | | X | | | | | | 0.048" hex |
| Dentsply Ankylos® C/X | X | | | | | | | | 1mm hex |
| Glidewell iNCLUSIVE® Tapered 3.7, 4.7, 5.2 | | | | | | | X | | 0.050"-1.25mm hex |
| Hiossen/Osstem HG mini (3.5) | | X | | | | | | | 1.2mm hex |
| Hiossen/Osstem HG standard (4.0, 4.5, 5.0) | | | | | X | | | | 1.2mm hex |
| Imtec® 3M Endure™ | | | | | | | X | | 0.050"-1.25mm hex |
| Implant Direct™ | | | | | X | | | | 0.050"-1.25mm hex |
| Keystone/Lifecore Genesis & Prima™ | | | | | X | | | | Square |
| Keystone/Lifecore Renova® | | | | | X | | | | 0.048" hex |
| Keystone/Lifecore Restore® | | | | | X | | | | 0.048" hex |
| Mega'Gen EZ Plus | | | | | | | X | | 0.048" hex |
| Mega'Gen Rescue | | | | | | | | X | 0.048" hex |
| Mega'Gen ExFeel internal (3.5) | | | | X | | | | | 1.2mm hex |
| Mega'Gen ExFeel internal (4.1, 4.8, 5.5) | | | | | | | X | | 1.2mm hex |
| Mega'Gen ExFeel external (3.3) | | | | X | | | | | 1.2mm hex |
| Mega'Gen ExFeel external (3.75, 4.0, 4.5, 5.0, 5.5) | | | | | | | X | | 1.2mm hex |
| MIS® Biocom & Seven | | | | | X | | | | 0.050"-1.25mm hex |
| Neoss | | | | | X | | | | Unigrip (not compatible with Nobel) |
| NobelActive™ & Replace™ | | | | | | | X | | Unigrip |
| NobelActive™ 3.0 | X | | | | | | | | Unigrip |
| OCO Biomedical | | | | X | | | | | 0.050"-1.25mm hex |
| Southern Implant External | | | | | | X | | | 0.048" hex |
| Southern Implant Tri Nex | | | | | | X | | | Unigrip |
| Southern Implant Octa | | | | | | X | | | Star/Torx |
| Straumann® Bone & Tissue Level | | | | | | | X | | Star/Torx |
| Sybron Pitt Easy | | | | X | | | | | 1.7mm hex |
| Thommen SPI® 3.5 | X | | | | | | | | Four lobe |
| Thommen SPI® 4.0, 4.5, 5.0, 6.0 | | | | X | | | | | Four lobe |
| Thommen SPI® zirconia abutment | | X | | | | | | | Four lobe |
| Zimmer® | | | | | X | | | | 0.050"-1.25mm hex |

NOTES

Atlantis™ titanium and zirconia abutments utilize the same torque and driver setting as the original implant manufacturer.
 Biomet 3i™ Certain® screws are gold plated only.
 Keystone Genesis and Prima™ abutment screws are made of titanium and have a titanium nitride coating.